STERILE EO

Pressure Products, Inc.
Customer Service
1861 N. Gaffey Street, Suite B
San Pedro, CA 90731 - USA
Tel +1-310-547-4973
Fax +1-310-547-4760
www.pressure-products.com

Oscor Inc.
3816 DeSoto Boulevard
Palm Harbor, FL 34683 - USA
Tel +1-727-937-2511
Fax +1-727-934-9835
www.oscor.com

EC REP


Osyka Medical GmbH
Albert-Einstein-Straße 3
12489 Berlin, Germany

For U.S.-California Only, Proposition 65, a State of California voter initiative, requires the following notice:
WARNING: This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Technical Manual

en-Instructions for Use

The SafeSept™ Transseptal Guidewire is a nitinol guidewire designed to cross the interatrial septum when supported by a transseptal introducer and be atraumatic when advanced unsupported into the left atrium.

A radiopaque coil is positioned on the distal end of the device to provide fluoroscopic localization during the procedure.

White printed markers are positioned on the proximal end of the device to provide approximate Transseptal Guidewire tip location relative to the needle tip during the procedure for use either with or without a hemostatic Y-adapter connected to the needle hub.

Intended Use

The SafeSept™ Transseptal Guidewire is used in conjunction with a transseptal needle to create the primary puncture in the interatrial septum and to guide the needle, dilator, and introducer through the septum from the right side of the heart to the left side.

Indications for use

Indicated for use in procedures where access to the left atrium via the transseptal technique is desired. The SafeSept™ Transseptal Guidewire is intended for single use only.

Contraindications

The use of the SafeSept™ Transseptal Guidewire is contraindicated in patients with the following conditions.

- Distorted anatomy due to congenital heart disease or other causes
- Significant chest or spine deformity
- The inability to lie flat
- Ongoing anticoagulation
- Left atrial thrombus or tumor
- Dilated aortic root
- Previous patch repair of the interatrial septum

Warnings and precautions

- For one use only, do not re-sterilize
- Store in a dry, dark, cool place
- Prolonged exposure to temperatures above 25°C (77°F) may damage the product
- Do not use if package is open or damaged in any way

Adverse events








In addition to all of the complications associated with any transseptal cardiac catheterization, the following can occur during the use of the SafeSept™ Transseptal Guidewire:

- Puncture of the atrial free wall
- Puncture of the aorta
- Puncture of the inferior vena cava
- Puncture of the coronary sinus
- Tamponade
- Hemothorax
- Arterial embolism from thrombus at the puncture site
- Stroke
- Death
- Atrial arrhythmias
- Residual atrial septal defects

Instructions for Use

Follow standard transseptal technique:

1. Per standard technique, aspirate and flush the transseptal system (including the needle, dilator, and sheath) before and after it is inserted into the right atrium.
2. Remove the protective tip cover and back the SafeSept™ Transseptal Guidewire slowly into the guidewire tip straightener. Do not pull the Transseptal Guidewire completely out of the tip straightener.
3. Fully insert the tip straightener into the hub of the transseptal needle.
4. Advance the SafeSept™ Transseptal Guidewire until the tip enters the curved portion of the transseptal needle (Printed markers are provided as a guide for common commercially available needle systems. See above. The Transseptal Guidewire can generally be advanced to the first printed marker when using a transseptal needle without a Y-adapter. Do not advance the guidewire beyond the tip of the needle).
5. Tent the fossa ovalis with the transseptal dilator. The transseptal needle tip position should be 0.5 to 1cm back from the tip of the dilator.
6. Under fluoroscopic guidance, **while tenting and maintaining constant force on the septum with the dilator**, slowly advance the SafeSept™ Transseptal Guidewire through the transseptal needle, dilator, and across the interatrial septum, into the left atrium. Continue to advance into one of the pulmonary veins (The radio-opaque coil on the guidewire should be seen to be within the left atrium and subsequently the pulmonary vein). **At no times should the guidewire be advanced when resistance is met without first determining the cause by fluoroscopy and taking remedial action. Confirm proper location of the SafeSept™ Transseptal Guidewire with fluoroscopy prior to proceeding.**
7. The Transseptal Guidewire should be advanced approximately 7cm into the left atrium so that the entire radiopaque coil is fluoroscopically visible to ensure that the transseptal needle tip does not engage the radiopaque coil during advancement.
8. Slowly advance the transseptal needle through the dilator and over the SafeSept™ Transseptal Guidewire across the septum. If desired, measure pressures and/or inject contrast through the needle before advancing the dilator and introducer sheath. Once the needle tip position is confirmed in the left atrium, advance the dilator and introducer sheath over the Transseptal Guidewire into the left atrium.
9. Slowly remove the SafeSept™ Transseptal Guidewire, needle, and dilator as a unit, leaving the introducer sheath in place. **Do not remove the Transseptal Guidewire first leaving the transseptal needle point exposed inside the left atrium. At no time should the introducer, needle, or guidewire be advanced or withdrawn when resistance is met without first determining the cause by fluoroscopy and taking remedial action.**
10. After use, the Transseptal Guidewire may be a potential biohazard. Handle and dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

en Instructions for Use	1
cs Návod k použití	2
da Brugsanvisning	2
de Gebrauchsanleitung	2
el Οδηγίες χρήσης	3
es Instrucciones para el uso	3
et Kasutusjuhend	4
fi Käyttöohjeet	4
fr Mode d'emploi	4
hu Használati utasítás	5
it Istruzioni per l'uso	5
ja 使用説明書	6
ko 사용 방법	6
lt Naudojimo instrukcijos	6
nl Gebruiksaanwijzing	7
no Bruksanvisning	7
pl Instrukcja użytkowania	8
pt Instruções de Utilização	8
ru Инструкции по использованию	9
sk Pokyny na používanie	9
sv Bruksanvisning	9
tr Kullanım Talimatları	10
zh 使用说明	10

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